



## Your reliable partner for aseptic filling

The filling of medical-pharmaceutical products is one of the most important of our sectors. As a **full-service provider** who contributes to people's well-being, we naturally set high quality and safety standards in filling, as well as in all other departments of our company.

FILLING



*Top quality manufacturing  
and filling of your syringes  
from one source.*



The aseptic filling, by which a sterile product is filled into vessels without the need for further terminal sterilisation, is a particularly critical process that not only requires qualified systems but also an extremely thorough validation of the aseptic conditions throughout the entire filling process. Transcoject is certified according to **ISO 13408** for aseptic filling and guarantees the permanent validation of all processes. We regularly review the quality and safety of the filling process by validating the procedure by means of integrated Media Fills that meet the current GMP requirements.

In the aseptic filling of medical devices, Transcoject has the experience and flexibility to meet all customer requirements. This way we realise the filling of low-volume lots up to batches of several hundred litres, depending on your **individual needs**. The production and filling of in-vitro medical devices is also part of the extensive range on offer at Transcoject.

Transcoject GmbH  
D-24539 Neumünster, Rügenstraße 8  
Fon +49 4321 885-0 • Fax +49 4321 885-65  
[www.transcoject.com](http://www.transcoject.com) • [info@transcoject.com](mailto:info@transcoject.com)